



# LactiLog™

Automated Data Collection System

**Sensational Solutions**  
for Today...and Tomorrow!

## LactiLog™ Automated Recording/Reporting System

A complete plug and play solution to automating data transfer and digitizing records, the LactiLog is custom-designed for unique reporting needs of the dairy industry. Simply connect the LactiCheck to your computer, download the LactiLog software and open the program – the results are then automatically recorded!

- ✓ Streamline data recording
- ✓ Easy access; transparent format
- ✓ Bar code compatible
- ✓ Eliminate costly transcription errors
- ✓ Instant availability of critical data facilitates management decisions



**Page & Pedersen International, Ltd.**

158 West Main Street Hopkinton, MA 01748 USA

[www.pagepedersen.com](http://www.pagepedersen.com)

# LactiLog™

## Automated Data Collection System

A streamlined & straight-forward solution to both data collection and reporting needs!



### Features and Benefits:

- ✓ Simple 'plug and play' installation. Connect the cable provided to your computer and the LactiCheck: Test recording begins automatically
- ✓ Readily integrated into existing systems
- ✓ Easy access, user configurable Suppliers and Products list for customized reporting needs
- ✓ Barcode compatible
- ✓ Tamper-proof data recording: Integrity of data is preserved through an encrypted file format for optimal security
- ✓ Sort by Supplier, Product or by Date
- ✓ Data is fully transferable thru comma separated (.csv) file format to standard spreadsheet programs for tracking and trending analysis, graphing and archiving
- ✓ Affordably introduces benefits of automated digitized data collection and real-time access to results for improving Process management analysis

Supplier	Date	Time	Fat	SNF	Density	Added Water	Freeze Point	Protein	Lactose	Product	Notes
LactiCal 3	5/24/2013	11:02:44 AM	3.46	9.27	31.48000	0.00	60.800	3.53	5.05	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:03:36 AM	3.47	9.29	31.55000	0.00	60.930	3.54	5.05	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:04:31 AM	3.46	9.30	31.60000	0.00	60.990	3.54	5.06	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:05:24 AM	3.45	9.28	31.53000	0.00	60.870	3.54	5.05	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:06:16 AM	3.46	9.30	31.60000	0.00	60.990	3.54	5.06	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:07:07 AM	3.46	9.29	31.56000	0.00	60.930	3.54	5.05	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:07:59 AM	3.46	9.30	31.60000	0.00	60.990	3.54	5.06	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:08:52 AM	3.46	9.29	31.56000	0.00	60.930	3.54	5.05	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:09:45 AM	3.46	9.30	31.60000	0.00	60.990	3.54	5.06	Whole Milk (Full Fat)	
LactiCal 3	5/24/2013	11:10:37 AM	3.45	9.28	31.53000	0.00	60.870	3.54	5.05	Whole Milk (Full Fat)	

### Specifications:

The LactiLog consists of a custom-designed software program in CD format and a dedicated RS-232 serial cable\*. The purpose of the System is to facilitate automated collection of test data generated by the LactiCheck™ Milk Analyzer, secure preservation of data integrity and provide easy-access to data for reporting requirements.

\* RS-232 to USB converter optional

### System Requirements:

**Hardware:** Pentium or higher, 20 MB free space on hard-drive, RS-232 output outlet\*, CD reading capabilities (built-in or auxiliary)

\*An RS-232 to USB converter can be provided

**Software:** Windows 98, 2000, ME, XP, Vista and Windows 7



## Page & Pedersen International, Ltd.

158 West Main Street, Hopkinton MA. 01748 U.S.A.

Tel: (508) 435 5966, Fax: (508) 435 8198

E-Mail: info@pagepedersen.com

[www.pagepedersen.com](http://www.pagepedersen.com)